

PRODUCT CATALOGUE

SSI-144 AUTOCLAVE (Portable)

Technical specifications

- Made of aluminum or stainless steel , having 63mm pressure guage , stop cock .
- Additional safety valve , steam release valve container and stand electrical.

Sizes available

Size (dia. X height)	Capacity (in liters)
300x300mm	20 Liters
300x375mm	25 Liters
300x450mm	30 Liters
300x500mm	40 Liters



SSI-145 AUTOCLAVE (Vertical) DOUBLE WALL

Technical specifications

- A versatile, electrically heated sterilizer under saturated steam of Adjustable pressure of 12 to 15 psi.
- Double walled construction with boiler & steam jacket made of thick stainless steel and gap is filled with high grade glass wool to bring the temperature loss to a minimal level.
- The outer chamber is made out of thick Stainless steel.
- The lid is made of thick stainless steel plate (joint less) and can be sealed by a joint less neoprene rubber gasket by wing nuts through hinged bolts by engaging slotted lugs on rim of the cover.
- Provided with wingnut (brass Nut) dial pressure gauge, air/ steam release cock, water level Indicator, spring loaded safety valve adjustable to any required point from 12 psi to 15 psi.
- A double safety valve has also been provided along with inlet.
- Supplied complete in Stainless Steel and Stainless steel dressing bins.

Sizes available

Size (dia. X height) in mm	Size (dia. X height) in inches	Capacity (in liters)	Load
250x450mm	10x18"	22 Liters	2 KW
300x500mm	12x20"	40 Liters	2 KW
350x550mm	14x22"	50 Liters	3 KW
400x600mm	16x24"	78 Liters	4 KW
450x600mm	18x24"	98 Liters	5 KW
550x750mm	22x30"	152 Liters	6 KW



S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C