

PRODUCT CATALOGUE

S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C

SSI-115 WATER BATH PRECISION (Constant Temp. control)

Technical specifications

- Double walled with inner chamber made of SS sheet (304 grade) and outer made of MS sheet duly epoxy powder coated.
- The gap between the walls is filled with special grade glass wool to avoid thermal losses.
- A glass window is provided on the front side for an easy inside view.
- A high speed stirrer with a motor of FHP is fitted to the unit for continuous stirring and thus maintaining a uniform temperature throughout.
- A drain plug is also fitted to the unit.
- The inner chamber fitted with immersion type heating element is covered with perforated removable diffuser tray.
- Temperature from ambient +5°C to 100°C, Accuracy $\pm 0.5^\circ\text{C}$ is controlled by Digital Temperature Indicator cum controller.
- Digital Timer and Alarm will also be provided.
- Supplied complete with cord and plug to work on 220/230 volts AC supply.



Sizes available

Inner Chamber Size (in mm)	No. of Holes	Capacity (in liters)
425x250x225mm	17x10x9"	24liters
450x375x375mm	18x15x15"	64liters
600x450x450mm	24x18x18"	122liters