

PRODUCT CATALOGUE

SSI-159 FUME HOOD

Construction

- Available in Wood, MS powder coated and Stainless Steel.
- Clean air Fume exhaust hoods used to drive-out hazardous chemical fumes, odors and gases generated during the laboratory process work.

Special Features

- Work surface made of polished granite, SS 304 or glazed ceramic tiles
- Counter balanced front sash with toughed glass view panels
- Available with wide range of accessories such as sink, water, inlet & outlet, gas/air, power points, optional FLP fittings, under bench cup-board etc.

Sizes available

Working Size (ft)	2' x 2' x 2'	3' x 2' x 2'	4' x 2' x 2'	5' x 2½' x 2½'
Model	SSI-159-2	SSI-159-3	SSI-159-4	SSI-159-5



S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C