

PRODUCT CATALOGUE

S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C

SSI-130 HEATING MANTLES

Technical specifications

- Body fabricated out of single piece aluminum and logged with special mineral wool .
- Provided with ON/OFF switch, pilot lamp, energy control knob etc.
- Provided with joint less hand knitted glass yarn capable with standing high temperature upto 400° C.
- Work on 220 VAC

Sizes available

Chamber Size	Wattage
50ml	60
100ml	60
250ml	150
500ml	200
1 liter	300
2 liter	450
3 liter	450
5 liter	450
10 liters	2x500
20 liters	4x600

