

PRODUCT CATALOGUE

SSI-131 HOT PLATE

Technical specifications

- Body fabricated out of single piece aluminum and lagged with special mineral wool .
- Provided with ON/OFF switch, pilot lamp, energy control knob etc.
- Provided with joint less hand knitted glass yarn capable with standing high temperature upto 400°C.
- Work on 220 VAC

Types of hot plate

- ROUND, SINGLE 20CM(8") DIA. Moulded top.
- ROUND, SINGLE 20CM(8") DIA. Moulded top.
- COIL TYPE, 20CM (8") DIA.



SSI-131D HOT PLATE (Rectangular)

Technical specifications

- Size : 300 x 450 mm (Approx.)
- Rectangular Hot Plate with Stainless Steel Top / Cast iron top
- Body: Made of thick mild steel sheet duly finished in powder coated paint
- Temperature Controlled by Energy Regulator
- Power : 220 Volts, Single Phase, 50 Hz.

Sizes available

Chamber Size	Wattage
30X25cm	12x10"
40x25cm	16x10"
45x30cm	18x12"
60x45cm	24x18"

