

PRODUCT CATALOGUE

SSI-157 LAMINAR AIR FLOW CABINET

Construction

- Available in SS 304, SS 316, MS, Melamine Laminated Wood materials.
- Horizontal / Vertical Laminar Airflow Benches allows operation in a sterile and particle free conditions. Due to continuous flushing of the working area by unidirectional and HEPA filtered airflow, the equipment assures a full product protection.

Special Features

International standard MINIPLEAT HEPA filters, high efficient & washable primary filters, Magnahelic gauges for accurate DP reading, spring suspended DIDW blower assembly for high performance at low sound level, SS working surface with curved front for convenient working.

Technical Specifications

Particle Retention	Air Velocity	Noise Level
0.3µ & above	90 FPM ± 20	65 DB ± 5

Sizes available

Working Size	2' x 2' x 2'	3' x 2' x 2'	4' x 2' x 2'	6' x 2' x 2'
Horizontal	SSI-157-H2	SSI-157-H3	SSI-157-H4	SSI-157-H6
Vertical	SSI-157-V2	SSI-157-V3	SSI-157-V4	SSI-157-V6



S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C