

PRODUCT CATALOGUE

SSI-111 SEROLOGICAL WATER BATH

Technical specifications

- Outer Body made of MS duly powder coated
- Inner Chamber made of SS 304
- Dimensions: L x W x D: 455 x 300 x 175 mm Capacity: 25 ltrs
- Temperature control from ambient to 110°C
- Temperature Accuracy: + 0.5°C
- With Microprocessor based PID controller
- Power: 220V, 50Hz.

Sizes available

Inner Chamber Size (in mm)	Inner Chamber Size (in inches)	Capacity (in liters)
300x175x175mm	12x7x7"	12liters
300x250x175mm	12x10x7"	14liters
375x300x175mm	15x12x7"	18liters
455x300x175mm	18x12x7"	25liters
600x300x175mm	24x12x7"	32liters



S
A
N
D
Y
S
C
I
E
N
T
I
F
I
C