

Ultra sonic-Size & Specification

VIBESONIC

Sl no	Model No	TUC2100	TUC3100	TUC5150	TUC6150
1	Capacity	2 Liters	3Liters	5 Liters (N)	6 Liters
2	Ultrasonic Frequency (Khz)	40	40	40	40
3	Ultrasonic Power	60W	100W	150W	200W
4	Ultrasonic Sensor	Piezo-PZT/USA	Piezo-PZT/USA	Piezo-PZT/USA	Piezo-PZT/USA
5	Mounting	Bottom Mounted	Bottom Mounted	Bottom Mounted	Bottom Mounted
6	Tank Internal Dimen(mm)Approx.	150*140*100	240*140*100	240*140*150	300*150*150
7	Heater With digital temperature control	Yes.100W	Yes.100W	Yes.150W	Yes.250W
8	Display	7 Segment	7 Segment	7 Segment	7 Segment
9	Digital Timer (0 to 99 min)	Yes	Yes	Yes	Yes
10	Degas Function	Yes	Yes	Yes	Yes
11	Drain Valve ½ ‘’	No	No	No	Yes
12	Tank Material	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
13	Outer Cabinet	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
14	Top Lid	Yes	Yes	Yes	Yes
15	Approx. Weight(kg)	4kg	4.8kg	6kg	9kg
16	Power Supply	230 VAC / 1 Phase with ground / Earth Connection			

VIBESONIC

Sl no	Model No	TUC10150	TUC14200	TUC20150	TUC30200
1	Capacity	9-10 Liters	14 Liters	20 Liters	30 Liters
2	Ultrasonic Frequency (Khz)	40	40	40	40
3	Ultrasonic Power	240W	300W	400W	600W
4	Ultrasonic Sensor	Piezo-PZT/USA	Piezo-PZT/USA	Piezo-PZT/USA	Piezo-PZT/USA
5	Mounting	Bottom Mounted	Bottom Mounted	Bottom Mounted	Bottom Mounted
6	Tank Internal Dimen(mm)Approx.	300*240*150	300*240*200	500*300*150	500*300*200
7	Heater With Digital Temperature control	Yes.400W	Yes.400W	Yes.600W	Yes.600W
8	Display	7 Segment	7 Segment	7 Segment	7 Segment
9	Digital Timer (0 to99 min)	Yes	Yes	Yes	Yes
10	Degas Function	Yes	Yes	Yes	Yes
11	Drain Valve ½ ‘	Yes	Yes	Yes	Yes
12	Tank Material	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
13	Outer Cabinet	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
14	Top Lid	Yes	Yes	Yes	Yes
15	Approx. Weight(kg)	11kg	12kg	18kg	20kg
16	Power Supply	230 VAC / 1 Phase with ground / Earth Connection			